



TALON[®]

LOW & MID-TOP PEDICLE SCREW

**FEATURES
AND
BENEFITS**

THORACOLUMBAR
FIXATION



OMNI SPINE
AN AMENDIA LINE

AMENDIA.

An all-Inclusive system, all for you.

The Talon® Low and Mid-Top Pedicle screw's comprehensive, multifaceted system easily changes and adjusts to the individual patient and surgeon's needs. Its top-loading, double lead thread design preserves spinal stability in skeletally mature patients as an adjunct to spinal fusion of the thoracic, lumbar, and/or sacral spine.

The Talon Pedicle Screw features:

- Improved preservation of spinal stability
- Double-lead thread screws for fast insertion and greater pull out strength
- 5.5mm system with different rod configurations
- 70 degree polyaxial design
- Ti-6Al-4V ELI alloy construction



TALON | LOW & MID-TOP PEDICLE SCREW

SOLID, LOW TOP-EXTENDED TAB

PART NUMBER	SIZE (D X L)
T04525	4.5 mm x 25 mm
T04530	4.5 mm x 30 mm
T04535	4.5 mm x 35 mm
T04540	4.5 mm x 40 mm
T04545	4.5 mm x 45 mm
T05030	5.0 mm x 30 mm
T05035	5.0 mm x 35 mm
T05040	5.0 mm x 40 mm
T05045	5.0 mm x 45 mm
T05050	5.0 mm x 50 mm
T05055	5.0 mm x 55 mm
T06030	6.0 mm x 30 mm
T06035	6.0 mm x 35 mm
T06040	6.0 mm x 40 mm
T06045	6.0 mm x 45 mm
T06050	6.0 mm x 50 mm
T06055	6.0 mm x 55 mm
T07030	7.0 mm x 30 mm
T07035	7.0 mm x 35 mm
T07040	7.0 mm x 40 mm
T07045	7.0 mm x 45 mm
T07050	7.0 mm x 50 mm
T07055	7.0 mm x 55 mm
T08030	8.0 mm x 30 mm
T08035	8.0 mm x 35 mm
T08040	8.0 mm x 40 mm
T08045	8.0 mm x 45 mm
T08050	8.0 mm x 50 mm
T08055	8.0 mm x 55 mm

SOLID, LOW-TOP

PART NUMBER	SIZE (D X L)
ST05035	5.0 mm x 35 mm
ST05040	5.0 mm x 40 mm
ST05045	5.0 mm x 45 mm
ST05050	5.0 mm x 50 mm
ST06035	6.0 mm x 35 mm
ST06040	6.0 mm x 40 mm
ST06045	6.0 mm x 45 mm
ST06050	6.0 mm x 50 mm
ST06055	6.0 mm x 55 mm
ST07035	7.0 mm x 35 mm
ST07040	7.0 mm x 40 mm
ST07045	7.0 mm x 45 mm
ST07050	7.0 mm x 50 mm
ST07055	7.0 mm x 55 mm
ST08035	8.0 mm x 35 mm
ST08040	8.0 mm x 40 mm
ST08045	8.0 mm x 45 mm
ST08050	8.0 mm x 50 mm

MID-TOP CANNULATED PEDICLE SCREW

PART NUMBER	SIZE (D X L)
CT05035	5.0 mm x 35 mm
CT05040	5.0 mm x 40 mm
CT05045	5.0 mm x 45 mm
CT05050	5.0 mm x 50 mm
CT05055	5.0 mm x 55 mm
CT06035	6.0 mm x 35 mm
CT06040	6.0 mm x 40 mm
CT06045	6.0 mm x 45 mm
CT06050	6.0 mm x 50 mm
CT06055	6.0 mm x 55 mm
CT07035	7.0 mm x 35 mm
CT07040	7.0 mm x 40 mm
CT07045	7.0 mm x 45 mm
CT07050	7.0 mm x 50 mm
CT07055	7.0 mm x 55 mm



For more information or to place an order call **1-877-755-3329**,
email customerservice@amendia.com or visit amendia.com.