

PIRANHHA

ANTERIOR CERVICAL PLATE

FEATURES AND BENEFITS

CERVICAL
FIXATION



OMNI SPINE
AN AMENDIA LINE

AMENDIA.

It's both thin and massive at the same time.

The Piranha Anterior Cervical Plate features a thin 2.2mm profile with large graft windows to enhance visualization and minimize tissue damage. The Piranha is anatomically designed with a small degree of lordosis and enhanced load sharing to promote fusion. The Piranha's unique locking mechanism provides tactile feedback and minimizes the possibility of screw loosening.

The Piranha features:

- Implant grade Titanium Alloy (Ti 6Al-4V ELI) construction
- Low-profile (2.2mm)
- Self-drilling Variable and Fixed screws in 12, 14, 16 and 18mm
- 16° cephalad/caudal screw angle with variable screws
- Large graft windows
- Lordosed plate design
- Available plates in 1, 2, 3 and 4 level sizes from 10-95mm



PIRANHA | ANTERIOR CERVICAL PLATE

1-LEVEL

| PART NUMBER | SIZE |
|-------------|------|
| P1-0010 | 10mm |
| P1-0012 | 12mm |
| P1-0014 | 14mm |
| P1-0016 | 16mm |
| P1-0018 | 18mm |
| P1-0020 | 20mm |
| P1-0022 | 22mm |
| P1-0024 | 24mm |

3-LEVEL

| PART NUMBER | SIZE |
|-------------|------|
| P3-0046 | 46mm |
| P3-0049 | 49mm |
| P3-0052 | 52mm |
| P3-0055 | 55mm |
| P3-0058 | 58mm |
| P3-0061 | 61mm |
| P3-0064 | 64mm |
| P3-0067 | 67mm |

2-LEVEL

| PART NUMBER | SIZE |
|-------------|------|
| P2-0024 | 24mm |
| P2-0026 | 26mm |
| P2-0028 | 28mm |
| P2-0030 | 30mm |
| P2-0032 | 32mm |
| P2-0035 | 35mm |
| P2-0038 | 38mm |
| P2-0041 | 41mm |
| P2-0044 | 44mm |
| P2-0047 | 47mm |

4-LEVEL

| PART NUMBER | SIZE |
|-------------|------|
| P4-0060 | 60mm |
| P4-0065 | 65mm |
| P4-0070 | 70mm |
| P4-0075 | 75mm |
| P4-0080 | 80mm |
| P4-0085 | 85mm |
| P4-0090 | 90mm |
| P4-0095 | 95mm |



For more information or to place an order call **1-877-755-3329**, email customerservice@amendia.com or visit amendia.com.